

**Human CLU (a.a.228–449)**

Swiss-Prot #	P10909
Protein names	CLU (APOJ) (CLI) (KUB1) (AAG4)
Gene names	
Organism	Human

Codon optimized cDNA Clone


cDNA Cat#	PEP-0728
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0728
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MMKTLLEFVGLLLTWESGQVLGDQTVSDNELQEMSNQGSKYVNKEIQNAVNGVKQIKTLIEKTNEERKTLNLEEAKK KKEDALNETRESETKLKLPGVCNETMMALWEECKPCLKQTCMKFYARVCRSGSLVGRQLEEFNQSSPFYFWMNGDR IDSLLENDRQQTHMLDVMQDHFSSRASSIIDELFQDRFFFTREPQDTYHYLPFSLPHRRPHFFFPKSRIVRSLMPFSPYEP LNFHAMFPFLEMIHEAQQAMDIFHSPAFQHPPTFEFIREGDDDRTVCREIRHNSTGCLRMKDQCDKREILSVDCSTN NPSQAKLRRELDSESLQVAERLTRKYNELLKSYQWKMLNTSSLEQLNEQFNWVSRLANLTQGEDQYYLRVTTVASHTSD SDVPSGVTEVVVKLFDSDPITVTVPVEVSRKNPKFMETVAEKALQEYRKKHREE

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**

**[Place Order>>>](#)**

Recombinant-Human-CLU-APOJ-CLI-KUB1-AAG4-Protein

recombinant protein, codon optimized, cDNA, P10909, Human, CLU, APOJ, CLI, KUB1, AAG4,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN