

**Human Bcl-2-like protein 1 (a.a.1-210)**

Swiss-Prot #	Q07817
Protein names	Bcl-2-like protein 1 (Bcl2-L-1) (Apoptosis regulator Bcl-X)
Gene names	BCL2L1 (BCL2L) (BCLX)
Organism	Human

Codon optimized cDNA Clone

cDNA Cat#	PEP-0492
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid


Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0492
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MSQSNRELVDVFLSYKLSQKGYSSQFSDVEENRTEAPEGTESEMETPSAINGNPSWHLADSPAVNGATGHSSSLDARE VIPMAAVKQALREAGDEFELRYRRAFSDLTSQLHITPGTAYQSFEQVVNELFRDGVNWGRIVAFFSFGGALCVESVDKE MQVLVSRIAAMMATYLNHLEPWIQENGWDTFVELYGNNAAAESRKGQERFNRWFLTGMTVAGVVLGSLFSRK

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**

**Place Order>>>**

Recombinant-Human-Bcl-2-Bcl2-L-cDNA-BCL2L1-BCL2L-BCLX-Bcl2-protein

recombinant protein, codon optimized, cDNA, Q07817, Human, BCL2L1, BCL2L, BCLX, Bcl-2-like protein, 1, Bcl2-L-1, Apoptosis regulator, Bcl-X,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN