

**Human Corticoliberin (a.a.25–153)**

Swiss-Prot #	P06850
Protein names	Corticoliberin (Corticotropin-releasing hormone) (Corticotropin-releasing factor) (CRF)
Gene names	CRH
Organism	Human

Codon optimized cDNA Clone

cDNA Cat# PEP-0457

Description Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat# PE-0457

Source Recombinant protein derived from E.coli. Or plant cell

Tag A 6x histidine tag

Formulation Liquid

Purity > 90%, as determined by SDS-PAGE under reducing conditions

Activity N/A


Applications Western Blot, ELISA or may be used for other applications determined by user

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Sequence MRLPLLVSAGVLLVALLPCPPCRALLSRGVPVGARQAPQHPQPLDFFQPPPPQSEQPQQPQARPVLLRMGEEYFLRLGNL
NKSPAAPLSPASSLLAGGSGSRPSPQATANFFRVLLQQLLLPRSLDSPAALAERGARNALGGHQEAPERERRSEPP
ISLDLTFHLLREVLEMARAEQLAQQAHSNRKLMEIIGK

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

Thousands
Recombinant
Protein, cDNA
Available>>>

**Place Order>>>**

Recombinant-Human-Corticoliberin-CRF-Protein-cDNA-CRH-CRF-Protein

recombinant, protein, codon, optimized, cDNA, P06850, Human, CRH, Corticoliberin, Corticotropin-releasing, hormone, Corticotropin-releasing, factor, CRF,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN