

**Human PI16 (a.a.28-442)**

|               |                                                  |
|---------------|--------------------------------------------------|
| Swiss-Prot #  | Q6UXB8                                           |
| Protein names | PI16 (CRISP9) (PSPBP) (PSEC0164) (UNQ289/PRO328) |
| Gene names    |                                                  |
| Organism      | Human                                            |

**Codon optimized cDNA Clone**

cDNA Cat# PEP-0743

**Description** Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

**Recombinant Protein**

Protein Cat# PE-0743

Source Recombinant protein derived from E.coli. Or plant cell

Tag A 6x histidine tag

Formulation Liquid

Purity &gt; 90%, as determined by SDS-PAGE under reducing conditions

Activity N/A


Applications Western Blot, ELISA or may be used for other applications determined by user

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

**Sequence** MHGSCSFLMLLPLLLLLLVATTGPVGALTDEEKRLMVELHNLRYAQVSP TASDMLHMRWDEELAAFAKAYARQCVWGHN  
KERGRGENLFAITDEGMDVPLAMEEWHHEREHYNLSAATCSPGQMGHYTQVVWAKTERIGCGSHFCEKLQGVETNI  
ELLVCNIEPPGNVKGKRPYQEGTPCSQCPSGYHCKNSLCEPIGSPEDAQDLPLYLVTEAPSFRAEASDSRKMGTTPSSLA  
TGIPAFVLVTEVSGSLATKALPAVETQAPTSLATKDPSPMATEAPPCVTTTEVPSILAAHSLPSLDEEPVTFPKSTHVPI  
KSADKVTDKTKVPSRSPENSLDPKMSLTGARELLPHAQEEAEAEALPPSSEVLASVFPQDKPGELQATLDHTGHTSS  
KSLPNFPNTSATANATGGRALALQSSLPGAEGPDKPSVVSGLNSGPGHVWGPLLGLLLLPLVLGIF

**Magnetic Beads**

Fast purifying  
DNA, RNA, Peptide  
protein, Antibody  
More>>>



Ion Exchange  
WCX, SCX, WAX  
Affinity Beads  
NH<sub>2</sub>, COOH, CHO, NHS,  
iodoacetyl, SH, OH,  
Silica, Thiol,  
etc. Conjugate Any  
Molecules C<sub>4</sub>, C<sub>8</sub>,  
C<sub>18</sub>, Endotoxin-  
removal, Bitin, Protein  
A, G, A/G, L,  
Streptavidin,  
Monoavidin GST  
More>>>

**Thousands  
Recombinant  
Protein, cDNA  
Available>>>**

**Place Order>>>**

Recombinant-Human-PI16-CRISP9-PSPBP-Protein-cDNA-PI16-CRISP9

recombinant protein, codon optimized, cDNA, Q6UXB8, Human, PI16, CRISP9, PSPBP, PSNQ289/PRO328,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN