

**Human ANO1 (a.a.628–731,813–882)**

Swiss-Prot #	Q5XXA6
Protein names	ANO1 (DOG1) (ORAOV2) (TAOS2) (TMEM16A)
Gene names	
Organism	Human
Codon optimized cDNA Clone	
cDNA Cat#	PEP-0469
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Recombinant Protein	
Protein Cat#	PE-0469
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	FKGRFVGRPGDYVYIFRSFRMEECAPGGCLMELCIQLSIIMLGKQLIQNNLFEIGIPKMKKLIRFKGRFVGRPGDYVYIFRSFRMEECAPGGCLMELCIQLSIIMLGKQLIQNNLFEIGIPKMKKLIRYLKLLKQQSPPDHEECVKKRQRYEVDYNLEPFAGLTPEYMEFAGLTPEYMEYLLKLLKQQSPPDHEECVKKRQRYEVDYNLEPFAGLTPEYME

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, Iodoacetyl, SH, OH, Silica, Thiol, etc. Conjugate Any Molecules C₄, C₈, C₁₈, Endotoxin-removal, Bitin, Protein A, G, A/G, L, Streptavidin, Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**



[Place Order>>>](#)

Recombinant-Human-ANO1-DOG1-ORAOV2-TAOS2-TMEM16A-Protein-cDNA

recombinant, protein, codon, optimized, cDNA, Q5XXA6, Human, ANO1, DOG1, ORAOV2, TAOS2, TMEM16A,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN