

**Human Glutathione S-transferase P (a.a.2-210)**

Swiss-Prot #	P09211
Protein names	Glutathione S-transferase P (EC 2.5.1.18) (GST class-pi) (GSTP1-1)
Gene names	GSTP1 (FAEES3) (GST3)
Organism	Human

Codon optimized cDNA Clone

cDNA Cat# PEP-0470

Description Applications: 1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat# PE-0470

Source Recombinant protein derived from E.coli. Or plant cell

Tag A 6x histidine tag

Formulation Liquid

Purity > 90%, as determined by SDS-PAGE under reducing conditions

Activity N/A


Applications Western Blot, ELISA or may be used for other applications determined by user

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Sequence MPPYTVVYFPVRGRCAALRMLLADQGQSWKEEVVTVETWQEGSLKASCLYQQLPKFQDGDLTLYQSNITLRHLGRTLGL
YGKDDQEEAALVDMVNDGVEDLRCKYISLIYTNYEAGKDDYVKALPGQLKPFETLLSQNGGKTFIVGDQISFADYNLLD
LLLIHEVLAPGCLDAFPLLSAYVGRLSARPKLKAFLASPEYVNLPIGNGKQ

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS,
iodoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**

**Place Order>>>**

Recombinant-Human-GST-GSTP1-1-Protein-cDNA-GSTP1-FAEES3-GST3

recombinant protein, codon optimized, cDNA, P09211, Human, GSTP1, FAEES3, GST3, Glutathione S-transferase P, GST, class-pi, GSTP1-1,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN