





**Human ANGPTL4 (ARP4) (HFARP) (PGAR) (PP1158) (PSEC0166)
(UNQ171/PRO197) (a.a. 26-406)**

Swiss-Prot #	Q9BY76
Protein names	Angiopoietin-related protein 4 (Angiopoietin-like 4) (Hepatic fibrinogen/angiopoietin-related protein) (HFARP)
Gene names	ANGPTL4 (ARP4) (HFARP) (PGAR) (PP1158) (PSEC0166) (UNQ171/PRO197)
Organism	Human
Codon optimized cDNA Clone	
cDNA Cat#	RPE-0169
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Recombinant Protein	
Protein Cat#	PE-0169
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MSGAPTAGAALMLCAATAVLLSAQGGPVQSKSPRFASWDEMNVLAHGLLQIQGLREHAE RTRSQLSALERRLSACGSACQGTGEGSTDLPLAPESRVDPEVLHSLQTQLKAQNSRIQQLF HKVAQQQRHLEKQHLRIQHLQSQFGLLDHKHLDHEVAKPARRKRLPEMAQPVDPAHNVSR LHRLPRDCQELFQVGERQSGLEIFQPGSPFFLVNCKMTSDGGWTVIQRHDSVDFNRP WEAYKAGFGDPHGEFVLGLEKVHSITGDRNSRLAVQLRDWDGNAELLQFSVHLGGEDTAY SLQLTAPVAGQLGATTVPVPSGLSVFSTWDQDHDLRDKNCAKSLSGGWVFGTCSHSNLT QGYFRSIPQQRQKLKKGIFWKTWRGRYYPLQATTMLIQPMAAEAAAS

Magnetic Beads



**Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>**



**Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Biotin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>**

**Thousands
Recombinant
Protein, cDNA
Available>>>**



Place Order >>>

Recombinant-protein-cDNA-Human-HFARP-ANGPTL4-ARP4-HFARP-PGAR-HFARP

recombinant, protein, codon, optimized, cDNA, Q9BY76, Human, ANGPTL4, ARP4, HFARP, PGAR, PP1158, PSNQ171/PRO197, Angiopoietin-related, protein, 4, Angiopoietin-like, 4, Hepatic,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN