

**Human TPBG (5T4) (a.a. 32-355)**

Swiss-Prot #	Q13641
Protein names	Trophoblast glycoprotein (5T4 oncofetal trophoblast glycoprotein) (5T4 oncotrophoblast glycoprotein) (5T4 oncofetal antigen) (M6P1)
Gene names	TPBG (5T4)
Organism	Human

Codon optimized cDNA Clone


cDNA Cat#	RPE-0182
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0182
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MPGGCSRGPAAAGDGRRLRLARLALVLLGWVSSSSPTSSASSFSSSAPFLASAVSAQPPLPD QCPALCECSEAARTVKCVNRNLTEVPTDLPAYVRNLFLTGNQLAVLPAGAFARRPPLAEL AALNLSGSRLEVRAGAFEHLPSLRQLDLSHNPLADLSPFAFSGSNASVSAPSPVLELIL NHIVPPEDERQNRSFEGMVVAALLAGRALQGLRRLLELASNHFLYLPDVLQALPSLRHLD LSNNSLVSLTYVSFRNLTHLESLHLEDNALKVLHNGTLAELQGLPHIRVFLDNNPWVDCD HMADMVTWLKETEVVQGKDRLTCAYPEKMRNRVLELNSADLDCDPILPPLSLQTSYVFLG IVLALIGAIFLLVLYLNRKGIKKWMHNIRDACRDHMEGYHYRYEINADPRLTNLSSNSDV

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
[More>>>](#)



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Biotin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
[More>>>](#)

Thousands
Recombinant
Protein, cDNA
Available>>>



[Place Order >>>](#)

Recombinant-protein-cDNA-Human-5T4-oncofetal-antigen-M6P1-TPBG-5T4

recombinant, protein, codon, optimized, cDNA, Q13641, Human, TPBG, 5T4, Trophoblast glycoprotein, 5T4, oncofetal, trophoblast glycoprotein, 5T4, oncotrophoblast glycoprotein, 5T4, oncofetal, antigen,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN