


**Human Alpha-galactosidase A (a.a.32-429)**

Swiss-Prot #	P06280
Protein names	Alpha-galactosidase A (EC 3.2.1.22) (Alpha-D-galactoside galactohydrolase) (Alpha-D-galactosidase A) (Melibiase) (Agalsidase)
Gene names	GLA
Organism	Human
Codon optimized cDNA Clone	
cDNA Cat#	PEP-0707
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Recombinant Protein	
Protein Cat#	PE-0707
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MQLRNPFLHIGCALALRFLALVSWDIPGARALDNLARTPTMGWLHWERFMCNLDCEPDSCISEKLFMEMAELMVSE GWKDAGYEYLCIDDCWMAPOQRDSEGLQADPQRFPHGIRQLANYVHSGKLGKIYADVGNKTCAGFPGSFGYYDIDAQT FADWGVDLLKFDGCYCDLENLADGYKHMSLALNRTGRSIVYSCEWPLYMWPFQKPNYTEIRQYCNHWRNFADIDDSWK SIKSILDWTSFNQERIVDVAGPGGWNDPDLVIGNFGLSWNQVTQMALWAIMAAPLFMSNDLRHISPOAKALLQDKDV IAINQDPLGKQGYQLRQGDNFEVWERPLSGLAWAVAMINRQEIGGPRSYTIAVASLGKGVACNPACFITQLLPVVRKLG FYEWTSRLRSHINPTGTVLLQLENTMQMSLKDLL

Magnetic Beads

**Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>**



**Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Biotin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>**

**Thousands
Recombinant
Protein, cDNA
Available>>>**



[Place Order>>>](#)

Recombinant-Human-Alpha-galactosidase-A-Melibiase-Protein-cDNA-GLA

recombinant protein, codon optimized, cDNA, P06280, Human, GLA, Alpha-galactosidase A, Alpha-D-galactoside galactohydrolase, Alpha-D-galactosidase A, Melibiase, Agalsidase,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN