


Human PARP3 (ADPRT3) (ADPRTL3) (a.a.313-533)

Swiss-Prot #	Q9Y6F1
Protein names	Poly [ADP-ribose] polymerase 3 (PARP-3) (EC 2.4.2.30) (NAD(+) ADP-ribosyltransferase 3) (Poly[ADP-ribose] synthetase 3) (pADPRT-3) (hPARP-3) (IRT1)
Gene names	PARP3 (ADPRT3) (ADPRTL3)
Organism	Human
Codon optimized cDNA Clone	
cDNA Cat#	RPE-0263
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.


Recombinant Protein

Protein Cat#	PE-0263
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MAPKPKPWVQTEGPEKKKGKQAGREEDPFRSTAEALKAIPA EKRI IRVDPTCPLSSNP GT QVYEDYNCTLNQTNIE NNNKFYIIQLLQDSNRFFTCWNHWGRVGEVGGQSKINHFRLED AKKDFEKKFKREKTKNNWAERDHFVSHPGKYTLIEVQAEDAEAEAVVKVDRGPVRTVTKRV QPCSLDPATQKLITNIFSKEMFKNTMALMDLDVKKMPLGKLSKQIARGFEALEALEEAL KGPTDGGQSLEELSSHFTYVIPHNFGHSQPPPI NSPELLQAKKDMLLVLADIELAQALQA VSEQKTVVEVPHPLDRDYQLLKCQLQLLD SGAP EYKVIQTYLEQTGSNHRCP TLQHIWK VNQEGEEDRFQAHSKLGNRKL LWHGTNMAVVAAILTSGLRIMPHSGGRVKGKIYFASENS KSAGYVIGMKCAHHVGYMFLGEVALGREHHINTDNPSLKSPPPGFDSVIARGHTEPDPT QDTELELDGQVVVPQGPVPCPEFSSSTFSQSEYLIYQESQCRLRYLLEVHL

Magnetic Beads



Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Biotin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

Thousands
Recombinant
Protein, cDNA
Available>>>



Place Order >>>

Recombinant-protein-cDNA-Human-PARP-3-pADPRT-hPARP-IRT1-PARP3

recombinant protein, codon optimized, cDNA, Q9Y6F1, Human, PARP3, ADPRT3, ADPRTL3, Poly[ADP-ribose] polymerase, 3, PARP-3, NAD(+) ADP-ribosyltransferase, 3, Poly[ADP-ribose] synthetase, 3, pADPRT-3, hPARP-3, IRT1

FOR IN VITRO RESEARCH USE ONLY. NOT FOR USE IN HUMAN