

**Human SLC45A3 (a.a.30-350)**

Swiss-Prot #	Q96JT2
Protein names	SLC45A3 (PCANAP6) (PRST)
Gene names	
Organism	Human

Codon optimized cDNA Clone

cDNA Cat#	PEP-0840
-----------	----------

Description Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0840
--------------	---------

Source Recombinant protein derived from E.coli. Or plant cell

Tag A 6x histidine tag

Formulation Liquid

Purity > 90%, as determined by SDS-PAGE under reducing conditions

Activity N/A


Applications Western Blot, ELISA or may be used for other applications determined by user

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Sequence MVQRLWVSRLRHRKAQLLLVNLTFGLFVCLAAAGITYVPPLLLEVGVVEEFMTMVLGIGPVLGLVCVPLLGASDHWGR
GRYGRRRPFIWALSIGILLSFLIPRAGWLAGLLCPDPRPLELALLILGVGLLDFCGQVCFTPLEALLSDLFRDPDHC
QAYSVYAFMISLGGCLGYLLPAIDWDTALAPYLGTEECFLGLTLIFLTCVAATLLVAEEAALGTEPAEGLSAPSL
SPHCCPCRARLAFRNLGALLPRLHQLCCRMPTLRRLFVAELCSWMALMTFTLFYTDVFGEGLYQGVPRAPGTEARRH
YDEGVRMGSLGLFLQCAISLVFSLVMDRLVQRFGTRAVYLASVAAPFVAAGATCLSHSVAVVTASAALTGFTFSALQIL
PYTLASLYHREKQVFLPKYRGDTGGASSEDLSMTSFLPGPKPGAPFPNGHVGAGGSGLLPPPPALCGASACDVSRRVV
GEPTEARVVPGRGICLDLAILDFAFLSQVAPSLFMGSIVQLSQSVTAYMVSAAGLGLVAIYFATQVVFDDKSLAKYSA

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**



[Place Order>>>](#)

Recombinant-Human-SLC45A3-PCANAP6-PRST-Protein-cDNA

recombinant protein, codon optimized, cDNA, Q96JT2, Human, SLC45A3, PCANAP6, PRST,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN