



Human PDGFD (IEGF) (SCDGFB) (MSTP036) (UNQ1899/PRO4345) (a.a.19-370)

Swiss-Prot #	Q9GZP0
Protein names	Platelet-derived growth factor D (PDGF-D) (Iris-expressed growth factor) (Spinal cord-derived growth factor B) (SCDGF-B) [Cleaved into: Platelet-derived growth factor D, latent form (PDGFD latent form); Platelet-derived growth factor D, receptor-binding f
Gene names	PDGFD (IEGF) (SCDGFB) (MSTP036) (UNQ1899/PRO4345)
Organism	Human

Codon optimized cDNA Clone

cDNA Cat# RPE-0148

Description Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat# PE-0148

Source Recombinant protein derived from E.coli. Or plant cell

Tag A 6x histidine tag

Formulation Liquid

Purity > 90%, as determined by SDS-PAGE under reducing conditions


Activity N/A

Applications Western Blot, ELISA or may be used for other applications determined by user


Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Sequence MHRLLFVYTLICANFCSCRDTSATPQSASIKALRNANLRDES NHLTDLYRRDETQVKG
NGYVQSPRFPNSYPRNLLLTWRLHSQENTRIQLVFDNQFGLLEAENDICRYDFVEVEDIS
ETSTIIRGRWCGHKEVPPRIKSRTNQIKITFKSDDYFVAKPGFKIYYSLEDFQPAAASE
TNWESVTSSISGVSYNSPVDPTLIADALDKKIAEFDTVEDLLKYFNPESWQEDLENMY
LDTPRYRGRSYHDKSKVDLRLNDDAKRYSTPRNYSVNIREELKLANVFFPRCLLVQ
RCGGNCGCGTVNWRSCNCGKTVKKYHEVLQFEPGHIKRRGRAKTALVDIQLDHHERC DCICSSRP

Magnetic Beads



Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Biotin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

Thousands
Recombinant
Protein, cDNA
Available>>>



Place Order >>>

Recombinant-protein-cDNA-Human-PDGF-D-PDGFD-IEGF-SCDGFB-MSTP036

recombinant protein, codon optimized, cDNA, Q9GZP0, Human, PDGFD, IEGF, SCDGFB, MSTP036, UNQ1899/PRO4345, Platelet-derived growth factor, D, PDGF-D, Iris-expressed growth factor, Spinal cord-

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN