

**Human CD84 (a.a.22-225)**


Swiss-Prot #	Q9UIB8
Protein names	CD84 (SLAMF5)
Gene names	
Organism	Human
Codon optimized cDNA Clone	
cDNA Cat#	PEP-0616
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0616
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MAQHHLWILLCLQTWPEAAGKDSEIFTVNGILGESVTFPVNIQEPRQVKIIAWTSKTSVAYVTPGDSETAPVVTVTHR NYYERIHAGPNYNLVISDLRMDAGDYKADINTQADPYTTTKRYNLQIYRRLGKPKITQSLMASVNSTCNVTLTCSVE KEEKNVTYNWSPLGEEGNVLQIFQTPEDQELTYTCTAQNPNVSNNSDSISARQLCADIAMGFRTHHTGLLSVLAMFFLLV LILSSVFLFRLFKRRQGRIFPEGSCLNTFTKNPYAASKKTIYTYIMASRNTQPAESRIYDEILQSKVLPSKEEPVNTVY SEVQFADKMGKASTQDSKPPGTSSYEIVI

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**

**Place Order>>>**

Recombinant-Human-CD84-SLAMF5-Protein-cDNA-CD84-SLAMF5-Protein

recombinant, protein, codon, optimized, cDNA, Q9UIB8, Human, CD84, SLAMF5,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN