

**Mouse Sparc (RP23-465I4.1-001) (a.a.17-302)**

Swiss-Prot #	Q5NCU5
Protein names	Putative uncharacterized protein (Secreted acidic cysteine rich glycoprotein) (Sparc protein)
Gene names	Sparc (RP23-465I4.1-001)
Organism	Mouse

Codon optimized cDNA Clone


cDNA Cat#	RPE-0310
Description	Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.
Tag	A 6x histidine tag
Formulation	Liquid
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0310
Source	Recombinant protein derived from E.coli. Or plant cell
Tag	A 6x histidine tag
Formulation	Liquid
Purity	> 90%, as determined by SDS-PAGE under reducing conditions
Activity	N/A
Applications	Western Blot, ELISA or may be used for other applications determined by user
Storage	-20° to -70° C. Avoid repeated freeze-thaw cycles.
Sequence	MRAWIFFLLCLAGRALAAPQQTEVAEEIVEEETVVEETGVPVGANPVQVEMGEFEDGAEETVEEVVADNCPQNHCKHGKVCDELDESNTPMCVCQDPTSCPAPIGEFEKVCSDNKTFFDSCHFFATKCTLEGTKKGKHLHDYIGPCKYIAPCLDSELTFFPLMRDWLKNVLTLYERDEGNNLLTEKQKLRLVKKIHENEKRLGADHPVELLARDFEKNYNMYIFPVHWQFGQLDQHPIDGYLSHTEPLAPLAPLIPMEHCTTRFFETCDLDNDKYIALEEWAGCFGIKEQDINKDL VI

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
[More>>>](#)



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, Iodoacetyl, SH, OH, Silica, Thiol, etc. Conjugate Any Molecules C4, C8, C18, Endotoxin-removal, Biotin, Protein A, G, A/G, L, Streptavidin, Monoavidin GST
[More>>>](#)

Thousands
Recombinant
Protein, cDNA
Available>>>



[Place Order >>>](#)

Recombinant-Sparc-cDNA-Sparc-Mouse-Sparc-protein-Sparc

recombinant protein, codon optimized, cDNA, Q5NCU5, Mouse, Sparc, RP23-465I4.1-001, Putative, uncharacterized, protein, Sacidic, cysteine, rich, glycoprotein, Sparc, protein,

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN