

**Human SPRY2 (a.a.1-315)**

Swiss-Prot #	O43597
Protein names	Protein sprouty homolog 2 (Spry-2)
Gene names	SPRY2
Organism	Human

Codon optimized cDNA Clone

cDNA Cat#	RPE-0101
-----------	----------

Description Applications:1. Codon-optimized cDNA sequence designed for recombinant protein production in both E.coli and Mammalian cell. The cDNA has been cloned into E.coli expression vector with 6x Histidine at N-terminus. 2. May be used as RNAi assistance to restore the loss-of-function.

Tag A 6x histidine tag

Formulation Liquid

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Recombinant Protein

Protein Cat#	PE-0101
--------------	---------

Source Recombinant protein derived from E.coli. Or plant cell

Tag A 6x histidine tag

Formulation Liquid

Purity > 90%, as determined by SDS-PAGE under reducing conditions

Activity N/A


Applications Western Blot, ELISA or may be used for other applications determined by user

Storage -20° to -70° C. Avoid repeated freeze-thaw cycles.

Sequence MEARAQSGNGSQPLLQTPRDGGRQGEPPDRDALTTQQVHVLSDQIRAIRNTNEYTEGPT
VVPRPGLKPAPRPSTQHKHERLHGLPEHRQPPRLQHSQVHSSARAPLSRSISTVSSGSR
STRSTSSSSSEQRLGSSFSFGPVADGIIRVQPKSELKPGELKPLSKEDLGLHAYRCED
CGKCKCKECTYPRPLPSDWICDKQCLCSAQNVIDYGTCCVCKGLFYHCSNDDNEDNCADN
PCSCSQSHCCTRWSAMGVMSLFLPLWCYLPAGKGLKLCQGCYDRVNRPGCRCKNSNTVC CKVPTVPPRNFEKPT

Magnetic Beads

Fast purifying
DNA, RNA, Peptide
protein, Antibody
More>>>



Ion Exchange
WCX, SCX, WAX
Affinity Beads
NH₂, COOH, CHO, NHS, I
odoacetyl, SH, OH,
Silica, Thiol,
etc. Conjugate Any
Molecules C₄, C₈,
C₁₈, Endotoxin-
removal, Bitin, Protein
A, G, A/G, L,
Streptavidin,
Monoavidin GST
More>>>

**Thousands
Recombinant
Protein, cDNA
Available>>>**

**Place Order >>>**

Recombinant-protein-cDNA-Human-Protein-sprouty-homolog-2-SPRY2

recombinant_protein_codon_optimized_cDNA_O43597_Human_SPRY2_Protein_sprouty_homolog_2_Spry-2

FOR in vitro RESEARCH USE ONLY. NOT FOR USE IN HUMAN